Флаг не определен Flag is not defined d  
multilang/variables  
images/medvision

**MedVision**

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---

### 4. Technical specifications

Name	Value
Size (length x width x height)	1300x500x400 mm
Weight	30 kg
Voltage	110 - 220 ( $\pm 10\%$ ) V
Single-phase frequency	50/60 Hz
Maximum power consumption	500 W
Accommodation area	4 m <sup>2</sup>
Minimum recommended work space size	3x3 m

## 4 TECHNICAL SPECIFICATIONS

---

Operating conditions	<p><b>Operation:</b></p> <p>Range of room temperature: 10 – 25 °C</p> <p>Range of room humidity: 30 – 75 % (non-condensing)</p> <p>Dust concentration in the air: Maximum single – 0.5 mg/m<sup>3</sup> Daily average – 0.15 mg/m<sup>3</sup></p> <p><b>Storage:</b></p> <p>Range of room temperature: 10 - 40 °C</p> <p>Range of room humidity: 30 - 75 % (non-condensing)</p> <p>Dust concentration in the air: Maximum single – 0.5 mg/m<sup>3</sup> Daily average – 0.15 mg/m<sup>3</sup></p>
Running time	4 hours
Break, minimum	30 min

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\_\_\_\_ / \_\_\_\_ / \_\_\_\_  
(Signature of the responsible person)

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---

## 8. Transportation

**Warning!** The simulator must not be subject to shock loads during handling operations.

### Land transportation

For land transportation, rigid package is used. The simulator is wrapped with air-bubble film and stretch film. Individual components (monitors, instruments, small parts) are packed in cardboard boxes. The packed equipment is mounted on a pallet with a wooden frame. This type of packaging provides necessary safety during land transportation and ensures cargo safety.

### Air transportation

For air transportation, original package is used. The simulator is wrapped with air-bubble film and stretch film. The wrapped equipment is packed into specially designed wood cases and locked. This type of packaging does not allow equipment to move in the cargo section of the aircraft and provides the necessary safety during transportation.

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## 10. Disposal & Recycling

### Packaging disposal and recycling

Packaging protects the simulator from damage during transportation. The packaging materials are environmentally friendly, easy to dispose of and recyclable. Recycling the package saves raw materials and reduces waste. We recommend delivering the packaging to a recycling center if possible.

### Simulator disposal and recycling

Electrical and electronic appliances often contain valuable materials. They also include certain substances, mixed components and parts necessary for the operation and safety of devices. If disposed of as household waste, or if handled improperly, they can cause damage to human health and the environment. So, never dispose of a worn out simulator as part of your household waste.

We recommend that you deliver the worn out simulator to a local entity accepting and recycling electrical and electronic devices. If any personal data is stored in the simulator subject to disposal and recycling, then it is your personal responsibility to delete it.

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---

## 12. Warranty

The warranty period is 12 months. The warranty period starts from the day of signing the Acceptance Report by the parties and can be extended by signing a supplementary contract with the manufacturer.

For the sake of clarity, keep documents confirming the transfer of goods and provide them to service staff when applying for warranty service.

During the warranty period, you can use your rights for free fixing of defects and repairs, provided that you use the product subject to its intended use and comply with the requirements for installation, connection and operation set forth in the manual and data sheet.

**Warning!** Each product has its unique serial number. The serial number fixes the date of manufacture and configuration of the product. When contacting the customer service, name the serial number of the product. This will help you to get faster and better service.

If, during the warranty period, the service department fixes the product malfunctions by replacing a part of it, the customer must return the removed part to the service staff at the end of the repair works (at the time the equipment malfunction is fixed). New parts, replaced during the repair process, may be of a color different from that of the original parts. Claims based on the replaced parts color would not be accepted.

The simulator hardware and software are being constantly improved. So, images used in the product data sheet and manuals may differ from the simulator current design. Claims based on images would not be accepted. We do our best to ensure the documentation

---

## 12 WARRANTY

---

contains the latest information about the simulator.

**The customer service is entitled to reject a warranty claim in the following cases:**

- Detected mechanical damage
  - Improper storage and/or careless transportation
  - Detected damage caused by inappropriate weather conditions during transportation, storage and operation
  - Detected effect of chemicals and moisture
  - Failure to comply with installation and connection procedures
  - Failure to comply with the user manual
  - Detected damage caused by failure to keep the equipment clean
  - Detected damage caused by improper use of consumables or use of care products not recommended by the manufacturer
  - Detected signs of disassembly, repair and other interventions by individuals not authorized to provide such services
  - Connecting to the power network of parameters not corresponding to the technical specifications of the simulator
  - Damage to the equipment caused by animals or insects
  - Illegal actions of the third parties (vandalism)
  - Force Majeure (fire, flood, natural disasters, etc.)
-

- Malfunctions caused by foreign objects, animals, insects and liquids found inside the equipment

**Warranty does not cover:**

- Works on the adjustment, cleaning and other care of the simulator specified in the user manual
- Replacement of consumables specified in the user manual.

**The following features are not defects:**

- Smell of new plastic or rubber observed during the first few days of operation
- Change in color and gloss of parts in the course of operation
- Noises (within the scope of sanitary standards) produced because of the way of operation of certain simulator components:
  - fan
  - water valve
  - electric relay
  - electric engine
  - compressor
- Noise caused by natural wear and tear (aging) of materials:
  - crackles during heating / cooling
  - squeaking

## 12 WARRANTY

---

- minor knocks of mobile mechanisms
- The need to replace consumables and wearing parts that have become unusable as a result of their natural wear and tear in the course of operation.

In case you need any service support, please find contacts in **Contact Information**

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© PROHP, Štefanec 10, 811 00 Bratislava